

INTRODUCING THE ALL NEW!

Unlimited cure depth, no microleakage,

bulKEZ PLUS™

does even more!



One step, no layering or capping saves time



Unlimited depth of cure



No additional primers or activators needed



Prevents gaps & leakage



Same great benefits as BulKEZ, PLUS:



Enhanced shade-matching, polishability and wear resistance



Do more! Crown fabrication, core build-ups, post cementation



Perfectly balanced flow, stackability and easy finishing

Introducing **bulKEZ PLUS™**

The game-changing dual-cure bulk fill composite just got even better, letting you do even more in less time!

Utilising the same IntellitTek™ Technology for unlimited cure depth and microleakage prevention as the original, BulKEZ PLUS now incorporates state of the art 80 nm spherical zirconia fillers, enhancing the aesthetics with better shade-matching and polishability and ensuring greater wear resistance. Handling has also been fine tuned to allow for a wider variety of cases, allowing you to do more in less time. BulKEZ PLUS offers results unlike any other composite from top to bottom!

Learn more about the bulk fill composite that's in a class of its own.
Freephone 0808 1000 888 or visit evident.co.uk




Evident
...making a difference

ZD ZEST DENTAL SOLUTIONS®

THE NEW AND IMPROVED **bulkeZ PLUS**™


- BulkeZ PLUS is an easy-to-place, dual-cure composite that combines unlimited cure depth with flowable cavity adaptation, high strength and enhanced durability—all in one, simple step.
- Self-cure IntelliTek™ Technology controls and directs shrinkage, eliminating leakage in all restorations.
- New 80nm spherical zirconia filler particles give market leading aesthetics, polishability and wear resistance.
- Handling and thixotropic properties have been precisely tuned to allow BulkeZ PLUS to be used in a wider variety of clinical situations.



Images courtesy of Dr. Troy Schmedding, Walnut Creek, CA


BulkeZ PLUS is in a Class of It's Own – Here are the Facts:

Conversion Map




BulkeZ PLUS reaches a higher degree of conversion than leading light-cured brands*

Gloss Wear & Resistance



BulkeZ PLUS offers greater glossiness, polish retention and wear resistance than leading brands*

BulkeZ PLUS vs Traditional Composite



Light-cured composites shrink towards the light, pulling away and creating leakage. BulkeZ PLUS eliminates microleakage.*

Ordering Information

PART NUMBER	DESCRIPTION
90-00171	BulkeZ PLUS A1/B1 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00172	BulkeZ PLUS A2/B2 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00173	BulkeZ PLUS A3/A3.5/B3 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00174	BulkeZ PLUS Bleach White Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00175	BulkeZ PLUS C2/C3 Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00176	BulkeZ PLUS Core White Refill, 1x 6 gm syringe, 6x tips (19 gauge)
90-00177	BulkeZ PLUS Trial Kit, 1x 6gm syringe (A1/B1), 1x 6gm syringe (A2/B2), 1x 6gm syringe (A3/A3.5/B3), 1x 6gm syringe (C2/C3), 1x 6gm syringe (BL/OP), 1x 6gm syringe (CORE White), 18x Tips (17 gauge), 18x Tips (19 gauge)
94792	BulkeZ Mixing Tips (17 Gauge), 25 pack
94793	BulkeZ Mixing Tips (19 Gauge), 25 pack

Evident is the sole distributor of BulkeZ PLUS in the UK and Ireland

Evident
...making a difference

Evident
Unit G10, Screenworks
22 Highbury Grove
London
N5 2EF

Freephone 0808 1000 888
sales@evident.co.uk
evident.co.uk